

Nail Cleaning Laser Device

Operating Instructions

1. Place the device on the nail and press the operating button.
2. Light will begin to blink and the device will emit visible blue light.
3. The device will operate for 7 minutes and shut off automatically.

Note: Device may not be used while charging battery. Please allow battery to be fully charged before use.

Usage instructions

It must be done on a daily basis.
It should be done at least once a day.



Safety

- ★ The device should be used in strict compliance with safety and operating precautions provided in this manual.
- ★ Please review the manual before initial use.
- ★ No side effects.
- ★ Direct or reflected laser radiation should be prevented from direct eye contact.
- ★ Pregnant women should not use the device close to womb.
- ★ Keep the device away from children.
- ★ Children under the age of 18 should be treated under the supervision of an adult.

Care and Cleaning

Exterior side of device may be cleaned with moistened cloth. Once in a few times of use it is recommended to extract the rubber padding located in the inner surface of the device and to clean it with alcohol. The device should not be handed to other people for use without prior cleaning of the rubber padding. Storage temperature: -20°C-50°C

Troubleshooting

Battery light is red	Rechargeable battery has run out of power	Connect the charger to charge it. It may not be operated while charging.
Device does not operate	Technical fault	Contact dealer/manufacturer

Technical Information

Laser wavelength: 905nm
Blue light wavelength: 470nm
Laser intensity: 25W
Maximum pulse intensity: 22W



Adaptor Specification

Voltage: 5V
Rating: 5VDC
Max: 0.5A
Charging time: 3-4 hours
Input: 100-240V 50/60Hz 0.2A